

## Product datasheet for **TP725692**

### UBE2Z Human Recombinant Protein

#### Product data:

Product Type:	Recombinant Proteins
Description:	Recombinant Human Ubiquitin-conjugating enzyme E2 Z(C-His)
Species:	Human
Expression cDNA Clone or AA Sequence:	T218-V354
Tag:	C-His
Buffer:	Supplied as a 0.2 um filtered solution of 500mM PBS, 300mM NaCl, 10%Glycerol, pH7.4
Note:	MADE TO ORDER   Recombinant Human Ubiquitin-conjugating enzyme E2 Z is produced by E.coli expression system and the target gene encoding T218-V354 is expressed with a His tag at the C-terminus.
Storage:	Store at $\leq -70^{\circ}\text{C}$ , stable for 6 months after receipt.Store at $\leq -70^{\circ}\text{C}$ , stable for 3 months under sterile conditions after opening.Please minimize freeze-thaw cycles.
Stability:	12 months from date of despatch
Locus ID:	65264
UniProt ID:	<a href="#">Q9H832</a>
Synonyms:	E2 ubiquitin-conjugating enzyme Z; Uba6-specific E2 conjugating enzyme 1; Use1; Ubiquitin carrier protein Z; Ubiquitin-protein ligase Z
Protein Pathways:	Acute myeloid leukemia, Pathways in cancer, Ubiquitin mediated proteolysis



[View online »](#)