

Product datasheet for **TP725693**

Myostatin Propeptide (MSTN) Mouse Recombinant Protein

Product data:

Product Type:	Recombinant Proteins
Description:	Recombinant Human GDF-8(C-His)
Species:	Mouse
Expression cDNA Clone or AA Sequence:	D101-Q320
Tag:	C-His
Buffer:	Supplied as a 0.2 um filtered solution of 500mM PBS, 300mM NaCl, 10%Glycerol, pH7.4
Note:	MADE TO ORDER Recombinant Human GDF-8 is produced by E.coli expression system and the target gene encoding D101-Q320 is expressed with a His tag at the C-terminus.
Storage:	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt.Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening.Please minimize freeze-thaw cycles.
Stability:	12 months from date of despatch
Locus ID:	2660
UniProt ID:	Q14793
Synonyms:	Growth/differentiation factor 8;Myostatin
Protein Families:	Druggable Genome, Secreted Protein



[View online »](#)