

## Product datasheet for **TP725697**

### Prdx3 Mouse Recombinant Protein

#### Product data:

|                                       |  |
|---------------------------------------|--|
| Product Type:                         | Recombinant Proteins   |
| Description:                          | Recombinant Human Prdx3 Protein(C-His)   |
| Species:                              | Mouse  |
| Expression cDNA Clone or AA Sequence: | A2-Q257  |
| Tag:                                  | C-His  |
| Buffer:                               | Supplied as a 0.2 um filtered solution of 500mM PBS, 300mM NaCl, 10%Glycerol, pH7.4  |
| Note:                                 | MADE TO ORDER   Recombinant Human Prdx3 Protein is produced by E.coli expression system and the target gene encoding A2-Q257 is expressed with a His tag at the C-terminus.                                  |
| Storage:                              | Store at $\leq -70^{\circ}\text{C}$ , stable for 6 months after receipt.Store at $\leq -70^{\circ}\text{C}$ , stable for 3 months under sterile conditions after opening.Please minimize freeze-thaw cycles. |
| Stability:                            | 12 months from date of despatch  |
| Locus ID:                             | 11757  |
| UniProt ID:                           | <a href="#">P20108</a>   |
| Synonyms:                             | Thioredoxin-dependent peroxide reductase, mitochondrial; Antioxidant protein 1; AOP-1;PRX III; Perioredoxin-3; Protein MER5; Thioredoxin-dependent peroxiredoxin 3   |



[View online »](#)