

Product datasheet for **TP725705**

Sonic Hedgehog (SHH) Mouse Recombinant Protein

Product data:

Product Type:	Recombinant Proteins
Description:	Recombinant Human Sonic hedgehog protein/SHH(K122R)(C-His)
Species:	Mouse
Expression cDNA Clone or AA Sequence:	C25-G198(K122R)
Tag:	C-His
Buffer:	Supplied as a 0.2 um filtered solution of 500mM PBS, 300mM NaCl, 10%Glycerol, pH7.4
Note:	MADE TO ORDER Recombinant Human Sonic hedgehog protein/SHH is produced by E.coli expression system and the target gene encoding C25-G198(K122R) is expressed with a His tag at the C-terminus.
Stability:	12 months from date of despatch
Locus ID:	6469
UniProt ID:	Q15465



[View online »](#)