

Product datasheet for **TP725707**

RUNX1T1 Human Recombinant Protein

Product data:

Product Type:	Recombinant Proteins
Description:	Recombinant Human Protein CBFA2T1(C-His)
Species:	Human
Expression cDNA Clone or AA Sequence:	R201-N420
Tag:	C-His
Buffer:	Supplied as a 0.2 um filtered solution of 500mM PBS, 300mM NaCl, 10%Glycerol, pH7.4
Note:	MADE TO ORDER Recombinant Human Protein CBFA2T1 is produced by E.coli expression system and the target gene encoding R201-N420 is expressed with a His tag at the C-terminus.
Storage:	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt.Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening.Please minimize freeze-thaw cycles.
Stability:	12 months from date of despatch
Locus ID:	862
UniProt ID:	Q06455
Synonyms:	Cyclin-D-related protein; Eight twenty one protein; Protein ETO; Protein MTG8; Zinc finger MYND domain-containing protein 2
Protein Families:	Transcription Factors
Protein Pathways:	Acute myeloid leukemia, Pathways in cancer



[View online »](#)