

Product datasheet for **TP725709**

ARNTL2 Human Recombinant Protein

Product data:

Product Type:	Recombinant Proteins
Description:	Recombinant Human ARNTL2 Protein(C-His)
Species:	Human
Expression cDNA Clone or AA Sequence:	A2-G350
Tag:	C-His
Buffer:	Supplied as a 0.2 um filtered solution of 500mM PBS, 300mM NaCl, 10%Glycerol, pH7.4
Note:	MADE TO ORDER Recombinant Human ARNTL2 Protein is produced by E.coli expression system and the target gene encoding A2-G350 is expressed with a His tag at the C-terminus.
Storage:	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt.Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening.Please minimize freeze-thaw cycles.
Stability:	12 months from date of despatch
Locus ID:	56938
UniProt ID:	Q8WYA1
Synonyms:	Aryl hydrocarbon receptor nuclear translocator-like protein 2; Basic-helix-loop-helix-PAS protein MOP9; Brain and muscle ARNT-like 2; CYCLE-like factor; CLIF; Class E basic helix-loop-helix protein 6; bHLHe6; Member of PAS protein 9; PAS domain-containing protein 9
Protein Families:	Druggable Genome, Transcription Factors



[View online »](#)