

Product datasheet for **TP725712**

RAP1A Human Recombinant Protein

Product data:

Product Type:	Recombinant Proteins
Description:	Recombinant Human Ras-related protein Rap-1A(C-His)
Species:	Human
Expression cDNA Clone or AA Sequence:	R2-L184
Tag:	C-His
Buffer:	Supplied as a 0.2 um filtered solution of 500mM PBS, 300mM NaCl, 10%Glycerol, pH7.4
Note:	MADE TO ORDER Recombinant Human Ras-related protein Rap-1A is produced by E.coli expression system and the target gene encoding R2-L184 is expressed with a His tag at the C-terminus.
Stability:	12 months from date of despatch
Locus ID:	5906
UniProt ID:	P62834



[View online »](#)