

## Product datasheet for **TP725717**

### 53BP1 (TP53BP1) Human Recombinant Protein

#### Product data:

Product Type:	Recombinant Proteins
Description:	Recombinant Human TP53-binding protein 1 (C-His)
Species:	Human
Expression cDNA Clone or AA Sequence:	E1624-Y1800
Tag:	C-His
Buffer:	Supplied as a 0.2 um filtered solution of 500mM PBS, 300mM NaCl, 10%Glycerol, pH7.4
Note:	MADE TO ORDER   Recombinant Human TP53-binding protein 1 is produced by E.coli expression system and the target gene encoding E1624-Y1800 is expressed with a His tag at the C-terminus.
Storage:	Store at $\leq -70^{\circ}\text{C}$ , stable for 6 months after receipt.Store at $\leq -70^{\circ}\text{C}$ , stable for 3 months under sterile conditions after opening.Please minimize freeze-thaw cycles.
Stability:	12 months from date of despatch
Locus ID:	7158
UniProt ID:	<a href="#">Q12888</a>
Synonyms:	53BP1; p53-binding protein 1;p53BP1
Protein Families:	Druggable Genome, Transcription Factors



[View online »](#)