

## Product datasheet for **TP725718**

### Acid Phosphatase (ACP1) Human Recombinant Protein

#### Product data:

Product Type:	Recombinant Proteins
Description:	Recombinant Human LMW-PTPase(C-His)
Species:	Human
Expression cDNA Clone or AA Sequence:	A2-E140
Tag:	C-His
Buffer:	Supplied as a 0.2 um filtered solution of 500mM PBS, 300mM NaCl, 10%Glycerol, pH7.4
Note:	MADE TO ORDER   Recombinant Human LMW-PTPase is produced by E.coli expression system and the target gene encoding A2-E140 is expressed with a His tag at the C-terminus.
Stability:	12 months from date of despatch
Locus ID:	52
UniProt ID:	<a href="#">P24666</a>



[View online »](#)