

## Product datasheet for **TP725720**

### **DYRK1B Human Recombinant Protein**

#### **Product data:**

Product Type:	Recombinant Proteins
Description:	Recombinant Human DYRK1B Protein(C-His)
Species:	Human
Expression cDNA Clone or AA Sequence:	R38-H237
Tag:	C-His
Buffer:	Supplied as a 0.2 um filtered solution of 500mM PBS, 300mM NaCl, 10%Glycerol, pH7.4
Note:	MADE TO ORDER   Recombinant Human DYRK1B Protein is produced by E.coli expression system and the target gene encoding R38-H237 is expressed with a His tag at the C-terminus.
Storage:	Store at $\leq -70^{\circ}\text{C}$ , stable for 6 months after receipt.Store at $\leq -70^{\circ}\text{C}$ , stable for 3 months under sterile conditions after opening.Please minimize freeze-thaw cycles.
Stability:	12 months from date of despatch
Locus ID:	9149
UniProt ID:	<u><a href="#">Q9Y463</a></u>

[View online »](#)