

Product datasheet for **TP725724**

SATB1 Human Recombinant Protein

Product data:

Product Type:	Recombinant Proteins
Description:	Recombinant Human SATB1 (N-His-FLAG)
Species:	Human
Expression cDNA Clone or AA Sequence:	Met1-Asp763
Tag:	N-6His-FLAG
Buffer:	Lyophilized from a 0.2 um filtered solution of PBS, pH 7.4
Note:	Recombinant Human SATB1 is produced by our mammalian expression system and the target gene encoding native Met1-Asp763 is expressed with a His and a FLAG tag at the N terminus.
Stability:	12 months from date of despatch
Locus ID:	6304
UniProt ID:	Q01826
Summary:	SATB1 is a key factor when it comes to integrating higher-order chromatin architecture with gene regulation. SATB1 is involved in MHC class-I gene organization and chromatin organization in thymocytes.



[View online »](#)