

Product datasheet for **TP725728**

mtTFA (TFAM) Human Recombinant Protein

Product data:

Product Type:	Recombinant Proteins
Description:	Recombinant Human mtTFA(C-His)
Species:	Human
Expression cDNA Clone or AA Sequence:	S43-C246
Tag:	C-His
Buffer:	Supplied as a 0.2 um filtered solution of 500mM PBS, 300mM NaCl, 10%Glycerol, pH7.4
Note:	MADE TO ORDER Recombinant Human mtTFA is produced by E.coli expression system and the target gene encoding S43-C246 is expressed with a His tag at the C-terminus.
Storage:	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt.Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening.Please minimize freeze-thaw cycles.
Stability:	12 months from date of despatch
Locus ID:	7019
UniProt ID:	Q00059
Synonyms:	Transcription factor A, mitochondrial; Mitochondrial transcription factor 1; MtTF1; Transcription factor 6; TCF-6; Transcription factor 6-like 2
Protein Families:	Druggable Genome, Transcription Factors
Protein Pathways:	Huntington's disease



[View online »](#)