

## Product datasheet for **TP725729**

### F box protein 38 (FBXO38) Human Recombinant Protein

#### Product data:

Product Type:	Recombinant Proteins
Description:	Recombinant Human F-box only protein 38/FBXO38 (C-His)
Species:	Human
Expression cDNA Clone or AA Sequence:	D520-L922
Tag:	C-His
Buffer:	Supplied as a 0.2 um filtered solution of 500mM PBS, 300mM NaCl, 10%Glycerol, pH7.4
Note:	MADE TO ORDER   Recombinant Human F-box only protein 38/FBXO38 is produced by E.coli expression system and the target gene encoding D520-L922 is expressed with a His tag at the C-terminus.
Stability:	12 months from date of despatch
Locus ID:	81545
UniProt ID:	<a href="#">Q6PII6</a>



[View online »](#)