

## Product datasheet for **TP725735**

### CRM1 (XPO1) Human Recombinant Protein

#### Product data:

Product Type:	Recombinant Proteins
Description:	Recombinant Human Exportin 1/XPO1 (C-His)
Species:	Human
Expression cDNA Clone or AA Sequence:	Met1-Asp1071
Tag:	C-6His
Buffer:	Supplied as a 0.2 um filtered solution of PBS, pH 7.4
Note:	Recombinant Human XPO1 is produced by our mammalian expression system and the target gene encoding native Met1-Asp1071 is expressed with a 6-His tag at the C terminus.
Stability:	12 months from date of despatch
Locus ID:	7514
UniProt ID:	<a href="#">O14980</a>
Summary:	XPO1 mediates the export of unspliced or not completely spliced RNAs out of the nucleus from different viruses. Those viruses include HIV-1, HTLV-1, and influenza A.



[View online »](#)