

Product datasheet for **TP725738**

CHFR Human Recombinant Protein

Product data:

Product Type:	Recombinant Proteins
Description:	Recombinant Human E3 ubiquitin-protein ligase CHFR(C-His)
Species:	Human
Expression cDNA Clone or AA Sequence:	Q14-G128
Tag:	C-His
Buffer:	Supplied as a 0.2 um filtered solution of 500mM PBS, 300mM NaCl, 10%Glycerol, pH7.4
Note:	MADE TO ORDER Recombinant Human E3 ubiquitin-protein ligase CHFR is produced by E.coli expression system and the target gene encoding Q14-G128 is expressed with a His tag at the C-terminus.
Storage:	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt.Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening.Please minimize freeze-thaw cycles.
Stability:	12 months from date of despatch
Locus ID:	55743
UniProt ID:	Q96EP1
Synonyms:	Checkpoint with forkhead and RING finger domains protein; RING finger protein 196; RING-type E3 ubiquitin transferase CHFR
Protein Families:	Druggable Genome



[View online »](#)