

## Product datasheet for **TP725739**

### TLE 1 (TLE1) Human Recombinant Protein

#### Product data:

Product Type:	Recombinant Proteins
Description:	Recombinant Human Transducin-like enhancer protein 1/TLE2(C-His)
Species:	Human
Expression cDNA Clone or AA Sequence:	H11-N200
Tag:	C-His
Buffer:	Supplied as a 0.2 um filtered solution of 500mM PBS, 300mM NaCl, 10%Glycerol, pH7.4
Note:	MADE TO ORDER   Recombinant Human Transducin-like enhancer protein 1/TLE2 is produced by E.coli expression system and the target gene encoding H11-N200 is expressed with a His tag at the C-terminus.
Storage:	Store at $\leq -70^{\circ}\text{C}$ , stable for 6 months after receipt.Store at $\leq -70^{\circ}\text{C}$ , stable for 3 months under sterile conditions after opening.Please minimize freeze-thaw cycles.
Stability:	12 months from date of despatch
Locus ID:	7088
UniProt ID:	<a href="#">Q04724</a>
Synonyms:	E(Sp1) homolog; Enhancer of split groucho-like protein 1; ESG1
Protein Families:	Transcription Factors



[View online »](#)