

Product datasheet for **TP725740**

TRPS1 Human Recombinant Protein

Product data:

Product Type:	Recombinant Proteins
Description:	Recombinant Human Zinc finger transcription factor Trps1(C-His)
Species:	Human
Expression cDNA Clone or AA Sequence:	L9-H358
Tag:	C-His
Buffer:	Supplied as a 0.2 um filtered solution of 500mM PBS, 300mM NaCl, 10%Glycerol, pH7.4
Note:	MADE TO ORDER Recombinant Human Zinc finger transcription factor Trps1 is produced by E.coli expression system and the target gene encoding L9-H358 is expressed with a His tag at the C-terminus.
Stability:	12 months from date of despatch
Locus ID:	7227
UniProt ID:	Q9UHF7



[View online »](#)