

Product datasheet for **TP725742**

NTH1 (NTHL1) Human Recombinant Protein

Product data:

Product Type:	Recombinant Proteins
Description:	Recombinant Human hNTH1(C-His)
Species:	Human
Expression cDNA Clone or AA Sequence:	G31-L312
Tag:	C-His
Buffer:	Supplied as a 0.2 um filtered solution of 500mM PBS, 300mM NaCl, 10%Glycerol, pH7.4
Note:	MADE TO ORDER Recombinant Human hNTH1 is produced by E.coli expression system and the target gene encoding G31-L312 is expressed with a His tag at the C-terminus.
Stability:	12 months from date of despatch
Locus ID:	4913
UniProt ID:	P78549



[View online »](#)