

## Product datasheet for **TP725745**

### ADNP Human Recombinant Protein

#### Product data:

|                                       |                                                                                                                                                                                                                     |
|---------------------------------------|---------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Product Type:                         | Recombinant Proteins                                                                                                                                                                                                |
| Description:                          | Recombinant Human Activity-dependent neuroprotector homeobox protein(C-His)                                                                                                                                         |
| Species:                              | Human                                                                                                                                                                                                               |
| Expression cDNA Clone or AA Sequence: | L941-A1102                                                                                                                                                                                                          |
| Tag:                                  | C-His                                                                                                                                                                                                               |
| Buffer:                               | Supplied as a 0.2 um filtered solution of 500mM PBS, 300mM NaCl, 10%Glycerol, pH7.4                                                                                                                                 |
| Note:                                 | MADE TO ORDER   Recombinant Human Activity-dependent neuroprotector homeobox protein is produced by E.coli expression system and the target gene encoding L941-A1102 is expressed with a His tag at the C-terminus. |
| Storage:                              | Store at $\leq -70^{\circ}\text{C}$ , stable for 6 months after receipt.Store at $\leq -70^{\circ}\text{C}$ , stable for 3 months under sterile conditions after opening.Please minimize freeze-thaw cycles.        |
| Stability:                            | 12 months from date of despatch                                                                                                                                                                                     |
| Locus ID:                             | 23394                                                                                                                                                                                                               |
| UniProt ID:                           | <a href="#">Q9H2P0</a>                                                                                                                                                                                              |
| Synonyms:                             | Activity-dependent neuroprotective protein; ADNP                                                                                                                                                                    |
| Protein Families:                     | Transcription Factors                                                                                                                                                                                               |



[View online »](#)