

Product datasheet for **TP725751**

KDM4C Human Recombinant Protein

Product data:

Product Type:	Recombinant Proteins
Description:	Recombinant Human Lysine-specific demethylase 4C(C-His)
Species:	Human
Expression cDNA Clone or AA Sequence:	K381-K588
Tag:	C-His
Buffer:	Supplied as a 0.2 um filtered solution of 500mM PBS, 300mM NaCl, 10%Glycerol, pH7.4
Note:	MADE TO ORDER Recombinant Human Lysine-specific demethylase 4C is produced by E.coli expression system and the target gene encoding K381-K588 is expressed with a His tag at the C-terminus.
Storage:	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt.Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening.Please minimize freeze-thaw cycles.
Stability:	12 months from date of despatch
Locus ID:	23081
UniProt ID:	Q9H3R0
Synonyms:	Gene amplified in squamous cell carcinoma 1 protein; GASC-1 protein; JmjC domain-containing histone demethylation protein 3C; Jumonji domain-containing protein 2C; [histone H3]-trimethyl-L-lysine(9) demethylase 4C
Protein Families:	Druggable Genome, Transcription Factors



[View online »](#)