

Product datasheet for **TP725757**

FAM3B Human Recombinant Protein

Product data:

Product Type:	Recombinant Proteins
Description:	Recombinant Human FAM3B (C-His)
Species:	Human
Expression cDNA Clone or AA Sequence:	Glu30-Ser235
Tag:	C-6His
Buffer:	Lyophilized from a 0.2 um filtered solution of PBS, pH 7.4
Note:	Recombinant Human FAM3B is produced by our mammalian expression system and the target gene encoding native Glu30-Ser235 is expressed with a 6-His tag at the C-terminus.
Stability:	12 months from date of despatch
Locus ID:	54097
UniProt ID:	P58499
Summary:	FAM3B is a secreted glycoprotein that is a part of the cytokine like-protein family. Due to alternative splicing and alternate signal peptide cleavage sites, multiple isoforms are known. This protein is expressed by pancreatic islet alpha and beta cells and can also induce beta cell apoptosis.



[View online »](#)