

Product datasheet for **TP725761**

FOXA3 Human Recombinant Protein

Product data:

Product Type:	Recombinant Proteins
Description:	Recombinant Human HNF-3G(C-His)
Species:	Human
Expression cDNA Clone or AA Sequence:	R180-G340
Tag:	C-His
Buffer:	Supplied as a 0.2 um filtered solution of 500mM PBS, 300mM NaCl, 10%Glycerol, pH7.4
Note:	MADE TO ORDER Recombinant Human HNF-3G is produced by E.coli expression system and the target gene encoding R180-G340 is expressed with a His tag at the C-terminus.
Stability:	12 months from date of despatch
Locus ID:	3171
UniProt ID:	P55318



[View online »](#)