

Product datasheet for **TP725765**

JUNB Human Recombinant Protein

Product data:

Product Type:	Recombinant Proteins
Description:	Recombinant Human Transcription factor JunB (C-His)
Species:	Human
Expression cDNA Clone or AA Sequence:	C2-A200
Tag:	C-His
Buffer:	Supplied as a 0.2 um filtered solution of 500mM PBS, 300mM NaCl, 10%Glycerol, pH7.4
Note:	MADE TO ORDER Recombinant Human Transcription factor JunB is produced by E.coli expression system and the target gene encoding C2-A200 is expressed with a His tag at the C-terminus.
Stability:	12 months from date of despatch
Locus ID:	3726
UniProt ID:	P17275



[View online »](#)