

Product datasheet for **TP725769**

SAFB Human Recombinant Protein

Product data:

Product Type:	Recombinant Proteins
Description:	Recombinant Human Scaffold attachment factor B1/SAF-B(C-His)
Species:	Human
Expression cDNA Clone or AA Sequence:	I477-E705
Tag:	C-His
Buffer:	Supplied as a 0.2 um filtered solution of 500mM PBS, 300mM NaCl, 10%Glycerol, pH7.4
Note:	MADE TO ORDER Recombinant Human Scaffold attachment factor B1/SAF-B is produced by E.coli expression system and the target gene encoding I477-E705 is expressed with a His tag at the C-terminus.
Storage:	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt.Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening.Please minimize freeze-thaw cycles.
Stability:	12 months from date of despatch
Locus ID:	6294
UniProt ID:	Q15424
Synonyms:	SAF-B1; HSP27 estrogen response element-TATA box-binding protein; HSP27 ERE-TATA-binding protein
Protein Families:	Druggable Genome, Transcription Factors



[View online »](#)