

Product datasheet for **TP725778**

NFE2L3 Human Recombinant Protein

Product data:

Product Type:	Recombinant Proteins
Description:	Recombinant Human NFE2-related factor 3(C-His)
Species:	Human
Expression cDNA Clone or AA Sequence:	S347-K694
Tag:	C-His
Buffer:	Supplied as a 0.2 um filtered solution of 500mM PBS, 300mM NaCl, 10%Glycerol, pH7.4
Note:	MADE TO ORDER Recombinant Human NFE2-related factor 3 is produced by E.coli expression system and the target gene encoding S347-K694 is expressed with a His tag at the C-terminus.
Storage:	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt.Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening.Please minimize freeze-thaw cycles.
Stability:	12 months from date of despatch
Locus ID:	9603
UniProt ID:	Q9Y4A8
Synonyms:	Nuclear factor erythroid 2-related factor 3; NF-E2-related factor 3; Nuclear factor, erythroid derived 2, like 3
Protein Families:	Transcription Factors, Transmembrane



[View online »](#)