

Product datasheet for **TP725781**

UBE2E1 Human Recombinant Protein

Product data:

Product Type:	Recombinant Proteins
Description:	Recombinant Human Ubiquitin-conjugating enzyme E2 E1(C-His)
Species:	Human
Expression cDNA Clone or AA Sequence:	S2-T193
Tag:	C-His
Buffer:	Supplied as a 0.2 um filtered solution of 500mM PBS, 300mM NaCl, 10%Glycerol, pH7.4
Note:	MADE TO ORDER Recombinant Human Ubiquitin-conjugating enzyme E2 E1 is produced by E.coli expression system and the target gene encoding S2-T193 is expressed with a His tag at the C-terminus.
Storage:	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt.Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening.Please minimize freeze-thaw cycles.
Stability:	12 months from date of despatch
Locus ID:	7324
UniProt ID:	P51965
Synonyms:	E2 ubiquitin-conjugating enzyme E1; E2 ubiquitin-conjugating enzyme E1; UbCH6; Ubiquitin carrier protein E1; Ubiquitin-protein ligase E1
Protein Pathways:	Ubiquitin mediated proteolysis



[View online »](#)