

Product datasheet for **TP725786**

PAG608 (ZMAT3) Human Recombinant Protein

Product data:

Product Type:	Recombinant Proteins
Description:	Recombinant Human Zinc finger matrin-type protein 3(C-His)
Species:	Human
Expression cDNA Clone or AA Sequence:	I2-V289
Tag:	C-His
Buffer:	Supplied as a 0.2 um filtered solution of 500mM PBS, 300mM NaCl, 10%Glycerol, pH7.4
Note:	MADE TO ORDER Recombinant Human Zinc finger matrin-type protein 3 is produced by E.coli expression system and the target gene encoding I2-V289 is expressed with a His tag at the C-terminus.
Stability:	12 months from date of despatch
Locus ID:	64393
UniProt ID:	Q9HA38



[View online »](#)