

Product datasheet for **TP725792**

MTGR1 (CBFA2T2) Human Recombinant Protein

Product data:

Product Type:	Recombinant Proteins
Description:	Recombinant Human Protein CBFA2T2(C-His)
Species:	Human
Expression cDNA Clone or AA Sequence:	P31-E239
Tag:	C-His
Buffer:	Supplied as a 0.2 um filtered solution of 500mM PBS, 300mM NaCl, 10%Glycerol, pH7.4
Note:	MADE TO ORDER Recombinant Human Protein CBFA2T2 is produced by E.coli expression system and the target gene encoding P31-E239 is expressed with a His tag at the C-terminus.
Storage:	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt.Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening.Please minimize freeze-thaw cycles.
Stability:	12 months from date of despatch
Locus ID:	9139
UniProt ID:	O43439
Synonyms:	ETO homologous on chromosome 20; MTG8-like protein; MTG8-related protein 1;Myeloid translocation-related protein 1; p85
Protein Families:	Transcription Factors



[View online »](#)