

Product datasheet for **TP725804**

BCL9L Human Recombinant Protein

Product data:

Product Type:	Recombinant Proteins
Description:	Recombinant Human BCL9-like protein(C-His)
Species:	Human
Expression cDNA Clone or AA Sequence:	H39-V206
Tag:	C-His
Buffer:	Supplied as a 0.2 um filtered solution of 500mM PBS, 300mM NaCl, 10%Glycerol, pH7.4
Note:	MADE TO ORDER Recombinant Human BCL9-like protein is produced by E.coli expression system and the target gene encoding H39-V206 is expressed with a His tag at the C-terminus.
Stability:	12 months from date of despatch
Locus ID:	283149
UniProt ID:	<u>Q86UU0</u>



[View online »](#)