

Product datasheet for **TP725813**

Antizyme inhibitor 1 (AZIN1) Human Recombinant Protein

Product data:

Product Type:	Recombinant Proteins
Description:	Recombinant Human Antizyme inhibitor 1/AZI/AZI1(C-His)
Species:	Human
Expression cDNA Clone or AA Sequence:	A149-A448
Tag:	C-His
Buffer:	Supplied as a 0.2 um filtered solution of 500mM PBS, 300mM NaCl, 10%Glycerol, pH7.4
Note:	MADE TO ORDER Recombinant Human Antizyme inhibitor 1/AZI/AZI1 is produced by E.coli expression system and the target gene encoding A149-A448 is expressed with a His tag at the C-terminus.
Stability:	12 months from date of despatch
Locus ID:	51582
UniProt ID:	Q14977



[View online »](#)