

Product datasheet for **TP725816**

MAD1 (MAD1L1) Human Recombinant Protein

Product data:

Product Type:	Recombinant Proteins
Description:	Recombinant Human Mitotic spindle assembly checkpoint protein MAD1(C-His)
Species:	Human
Expression cDNA Clone or AA Sequence:	S597-A718
Tag:	C-His
Buffer:	Supplied as a 0.2 um filtered solution of 500mM PBS, 300mM NaCl, 10%Glycerol, pH7.4
Note:	MADE TO ORDER Recombinant Human Mitotic spindle assembly checkpoint protein MAD1 is produced by E.coli expression system and the target gene encoding S597-A718 is expressed with a His tag at the C-terminus.
Stability:	12 months from date of despatch
Locus ID:	8379
UniProt ID:	Q9Y6D9



[View online »](#)