

Product datasheet for **TP725823**

Human Recombinant Protein

Product data:

Product Type:	Recombinant Proteins
Description:	Recombinant Human Interferon alpha-2/IFNA2 (C-Fc)
Species:	Human
Expression cDNA Clone or AA Sequence:	Cys24-Glu188
Tag:	C-Fc
Buffer:	Lyophilized from a 0.2 um filtered solution of PBS, pH7.4
Note:	Recombinant Human IFNA2 is produced by our Mammalian expression system and the target gene encoding Cys24-Glu188 is expressed with a Fc tag at the C-terminus.
Stability:	12 months from date of despatch
Summary:	IFN-a is a part of the cytokine family and has a immunomodulatory function. This protein can activate dendritic cells and can be used for antiviral and anticancer treatment.



[View online »](#)