

Product datasheet for **TP725826**

SYCE1 Human Recombinant Protein

Product data:

Product Type:	Recombinant Proteins
Description:	Recombinant Human Synaptonemal complex central element protein 1(C-His)
Species:	Human
Expression cDNA Clone or AA Sequence:	E132-F351
Tag:	C-His
Buffer:	Supplied as a 0.2 um filtered solution of 500mM PBS, 300mM NaCl, 10%Glycerol, pH7.4
Note:	MADE TO ORDER Recombinant Human Synaptonemal complex central element protein 1 is produced by E.coli expression system and the target gene encoding E132-F351 is expressed with a His tag at the C-terminus.
Storage:	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt.Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening.Please minimize freeze-thaw cycles.
Stability:	12 months from date of despatch
Locus ID:	93426
UniProt ID:	<u>Q8N0S2</u>
Synonyms:	Cancer/testis antigen 76; CT76



[View online »](#)