

Product datasheet for **TP725835**

ASF1B Human Recombinant Protein

Product data:

Product Type:	Recombinant Proteins
Description:	Recombinant Human Histone chaperone ASF1B(C-His)
Species:	Human
Expression cDNA Clone or AA Sequence:	A2-I202
Tag:	C-His
Buffer:	Supplied as a 0.2 um filtered solution of 500mM PBS, 300mM NaCl, 10%Glycerol, pH7.4
Note:	MADE TO ORDER Recombinant Human Histone chaperone ASF1B is produced by E.coli expression system and the target gene encoding A2-I202 is expressed with a His tag at the C-terminus.
Storage:	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt.Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening.Please minimize freeze-thaw cycles.
Stability:	12 months from date of despatch
Locus ID:	55723
UniProt ID:	Q9NVP2
Synonyms:	Anti-silencing function protein 1 homolog B; hAsf1; hAsf1b; CCG1-interacting factor A-II; CIA-II; hCIA-II



[View online »](#)