

Product datasheet for **TP725837**

BRAF35 (HMG20B) Human Recombinant Protein

Product data:

Product Type:	Recombinant Proteins
Description:	Recombinant Human SMARCE1-related protein(C-His)
Species:	Human
Expression cDNA Clone or AA Sequence:	S2-L317
Tag:	C-His
Buffer:	Supplied as a 0.2 um filtered solution of 500mM PBS, 300mM NaCl, 10%Glycerol, pH7.4
Note:	MADE TO ORDER Recombinant Human SMARCE1-related protein is produced by E.coli expression system and the target gene encoding S2-L317 is expressed with a His tag at the C-terminus.
Storage:	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt.Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening.Please minimize freeze-thaw cycles.
Stability:	12 months from date of despatch
Locus ID:	10362
UniProt ID:	Q9P0W2
Synonyms:	SWI/SNF-related matrix-associated actin-dependent regulator of chromatin subfamily E member 1-related; BRCA2-associated factor 35; HMG box-containing protein 20B; HMG domain-containing protein 2; HMG domain-containing protein HMGX2; Sox-like transcriptional factor; Structural DNA-binding protein BRAF35
Protein Families:	Druggable Genome, Transcription Factors



[View online »](#)