

## Product datasheet for **TP725844**

### SIGLEC6 Human Recombinant Protein

#### Product data:

Product Type:	Recombinant Proteins
Description:	Recombinant Human SIGLEC-6 (C-His)
Species:	Human
Expression cDNA Clone or AA Sequence:	Gln27-Val331
Tag:	C-6His
Buffer:	Lyophilized from a 0.2 um filtered solution of PBS, pH7.4
Note:	Recombinant Human SIGLEC-6 is produced by our Mammalian expression system and the target gene encoding Gln27-Val331 is expressed with a His tag at the C-terminis
Stability:	12 months from date of despatch
RefSeq:	<a href="#">NP_942142</a>
Locus ID:	946
UniProt ID:	<a href="#">O43699</a>
Summary:	Siglec-6 is a sialic acid binding Ig-like lectin and I-type (Ig-type) lectin that belong to the Ig superfamily. This protein mediates sialic-acid binding to cells. This protein localizes in multiple compartments such as extracellular region and membrane function. This protein expression has been specualted to help slow the tempo of labor during childbirth. On the otherhand this protein expression upregulates pre-eclampsia.



[View online »](#)