

Product datasheet for **TP725850**

DNA polymerase delta p50 (POLD2) Human Recombinant Protein

Product data:

Product Type:	Recombinant Proteins
Description:	Recombinant Human DNA polymerase delta subunit 2/POLD2(C-His)
Species:	Human
Expression cDNA Clone or AA Sequence:	L85-L325
Tag:	C-His
Buffer:	Supplied as a 0.2 um filtered solution of 500mM PBS, 300mM NaCl, 10%Glycerol, pH7.4
Note:	MADE TO ORDER Recombinant Human DNA polymerase delta subunit 2/POLD2 is produced by E.coli expression system and the target gene encoding L85-L325 is expressed with a His tag at the C-terminus.
Stability:	12 months from date of despatch
Locus ID:	5425
UniProt ID:	P49005



[View online »](#)