

## Product datasheet for **TP725852**

### HEMK2 (N6AMT1) Human Recombinant Protein

#### Product data:

Product Type:	Recombinant Proteins
Description:	Recombinant Human Methyltransferase N6AMT1(C-His)
Species:	Human
Expression cDNA Clone or AA Sequence:	A2-S214
Tag:	C-His
Buffer:	Supplied as a 0.2 um filtered solution of 500mM PBS, 300mM NaCl, 10%Glycerol, pH7.4
Note:	MADE TO ORDER   Recombinant Human Methyltransferase N6AMT1 is produced by E.coli expression system and the target gene encoding A2-S214 is expressed with a His tag at the C-terminus.
Storage:	Store at $\leq -70^{\circ}\text{C}$ , stable for 6 months after receipt.Store at $\leq -70^{\circ}\text{C}$ , stable for 3 months under sterile conditions after opening.Please minimize freeze-thaw cycles.
Stability:	12 months from date of despatch
Locus ID:	29104
UniProt ID:	<a href="#">Q9Y5N5</a>
Synonyms:	HemK methyltransferase family member 2; M.HsaHemK2P; Lysine N-methyltransferase 9; Methylarsonite methyltransferase N6AMT1; Protein N(5)-glutamine methyltransferase
Protein Families:	Druggable Genome



[View online »](#)