

## Product datasheet for **TP725856**

### SF3A3 Human Recombinant Protein

#### Product data:

Product Type:	Recombinant Proteins
Description:	Recombinant Human Splicing factor 3A subunit 3/SF3A3(C-His)
Species:	Human
Expression cDNA Clone or AA Sequence:	E222-L501
Tag:	C-His
Buffer:	Supplied as a 0.2 um filtered solution of 500mM PBS, 300mM NaCl, 10%Glycerol, pH7.4
Note:	MADE TO ORDER   Recombinant Human Splicing factor 3A subunit 3/SF3A3 is produced by E.coli expression system and the target gene encoding E222-L501 is expressed with a His tag at the C-terminus.
Storage:	Store at $\leq -70^{\circ}\text{C}$ , stable for 6 months after receipt.Store at $\leq -70^{\circ}\text{C}$ , stable for 3 months under sterile conditions after opening.Please minimize freeze-thaw cycles.
Stability:	12 months from date of despatch
Locus ID:	10946
UniProt ID:	<a href="#">Q12874</a>
Synonyms:	SF3a60; Spliceosome-associated protein 61; SAP 61
Protein Pathways:	Spliceosome



[View online »](#)