

## Product datasheet for **TP725862**

### WDR4 Human Recombinant Protein

#### Product data:

Product Type:	Recombinant Proteins
Description:	Recombinant Human WD repeat-containing protein 4(C-His)
Species:	Human
Expression cDNA Clone or AA Sequence:	A28-L140
Tag:	C-His
Buffer:	Supplied as a 0.2 um filtered solution of 500mM PBS, 300mM NaCl, 10%Glycerol, pH7.4
Note:	MADE TO ORDER   Recombinant Human WD repeat-containing protein 4 is produced by E.coli expression system and the target gene encoding A28-L140 is expressed with a His tag at the C-terminus.
Storage:	Store at $\leq -70^{\circ}\text{C}$ , stable for 6 months after receipt.Store at $\leq -70^{\circ}\text{C}$ , stable for 3 months under sterile conditions after opening.Please minimize freeze-thaw cycles.
Stability:	12 months from date of despatch
Locus ID:	10785
UniProt ID:	<a href="#">P57081</a>
Synonyms:	tRNA (guanine-N(7)-)-methyltransferase non-catalytic subunit WDR4; Protein Wuho homolog; hWH;WD repeat-containing protein 4



[View online »](#)