

Product datasheet for **TP725866**

DACH2 Human Recombinant Protein

Product data:

Product Type:	Recombinant Proteins
Description:	Recombinant Human Dachshund homolog 2(C-His)
Species:	Human
Expression cDNA Clone or AA Sequence:	Q419-Q583
Tag:	C-His
Buffer:	Supplied as a 0.2 um filtered solution of 500mM PBS, 300mM NaCl, 10%Glycerol, pH7.4
Note:	MADE TO ORDER Recombinant Human Dachshund homolog 2 is produced by E.coli expression system and the target gene encoding Q419-Q583 is expressed with a His tag at the C-terminus.
Stability:	12 months from date of despatch
Locus ID:	117154
UniProt ID:	<u>Q96NX9</u>



[View online »](#)