

Product datasheet for **TP725871**

ATG4A Human Recombinant Protein

Product data:

Product Type:	Recombinant Proteins
Description:	Recombinant Human Cysteine protease ATG4A (C-His)
Species:	Human
Expression cDNA Clone or AA Sequence:	T289-V398
Tag:	C-His
Buffer:	Supplied as a 0.2 um filtered solution of 500mM PBS, 300mM NaCl, 10%Glycerol, pH7.4
Note:	MADE TO ORDER Recombinant Human Cysteine protease ATG4A is produced by E.coli expression system and the target gene encoding T289-V398 is expressed with a His tag at the C-terminus.
Stability:	12 months from date of despatch
Locus ID:	115201
UniProt ID:	Q8WYN0



[View online »](#)