

Product datasheet for **TP725875**

Angiopoietin like 7 (ANGPTL7) Human Recombinant Protein

Product data:

Product Type:	Recombinant Proteins
Description:	Recombinant Human Angiopoietin-related protein 7(C-His)
Species:	Human
Expression cDNA Clone or AA Sequence:	K28-G170
Tag:	C-His
Buffer:	Supplied as a 0.2 um filtered solution of 500mM PBS, 300mM NaCl, 10%Glycerol, pH7.4
Note:	MADE TO ORDER Recombinant Human Angiopoietin-related protein 7 is produced by E.coli expression system and the target gene encoding K28-G170 is expressed with a His tag at the C-terminus.
Storage:	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt.Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening.Please minimize freeze-thaw cycles.
Stability:	12 months from date of despatch
Locus ID:	10218
UniProt ID:	O43827
Synonyms:	Angiopoietin-like factor; Angiopoietin-like protein 7; Cornea-derived transcript 6 protein
Protein Families:	Druggable Genome, Secreted Protein, Transmembrane



[View online »](#)