

Product datasheet for **TP725877**

Thimet Oligopeptidase (THOP1) Human Recombinant Protein

Product data:

Product Type:	Recombinant Proteins
Description:	Recombinant Human Thimet oligopeptidase /THOP1(C-His)
Species:	Human
Expression cDNA Clone or AA Sequence:	P90-K220
Tag:	C-His
Buffer:	Supplied as a 0.2 um filtered solution of 500mM PBS, 300mM NaCl, 10%Glycerol, pH7.4
Note:	MADE TO ORDER Recombinant Human Thimet oligopeptidase /THOP1 is produced by E.coli expression system and the target gene encoding P90-K220 is expressed with a His tag at the C-terminus.
Stability:	12 months from date of despatch
Locus ID:	7064
UniProt ID:	P52888



[View online »](#)