

Product datasheet for **TP725880**

CKAP1 (TBCB) Human Recombinant Protein

Product data:

Product Type:	Recombinant Proteins
Description:	Recombinant Human Tubulin-folding cofactor B/TBCB(C-His)
Species:	Human
Expression cDNA Clone or AA Sequence:	E2-I244
Tag:	C-His
Buffer:	Supplied as a 0.2 um filtered solution of 500mM PBS, 300mM NaCl, 10%Glycerol, pH7.4
Note:	MADE TO ORDER Recombinant Human Tubulin-folding cofactor B/TBCB is produced by E.coli expression system and the target gene encoding E2-I244 is expressed with a His tag at the C-terminus.
Storage:	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt.Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening.Please minimize freeze-thaw cycles.
Stability:	12 months from date of despatch
Locus ID:	1155
UniProt ID:	Q99426
Synonyms:	Cytoskeleton-associated protein 1; Cytoskeleton-associated protein CKAP1; Tubulin-specific chaperone B



[View online »](#)