

Product datasheet for **TP725898**

VTI1A Human Recombinant Protein

Product data:

Product Type:	Recombinant Proteins
Description:	Recombinant Human VTI1A Protein(C-His)
Species:	Human
Expression cDNA Clone or AA Sequence:	S2-A203
Tag:	C-His
Buffer:	Supplied as a 0.2 um filtered solution of 500mM PBS, 300mM NaCl, 10%Glycerol, pH7.4
Note:	MADE TO ORDER Recombinant Human VTI1A Protein is produced by E.coli expression system and the target gene encoding S2-A203 is expressed with a His tag at the C-terminus.
Storage:	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt.Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening.Please minimize freeze-thaw cycles.
Stability:	12 months from date of despatch
Locus ID:	143187
UniProt ID:	Q96AJ9
Synonyms:	Vesicle transport through interaction with t-SNAREs homolog 1A; Vesicle transport v-SNARE protein Vti1-like 2; Vti1-rp2
Protein Families:	Transmembrane
Protein Pathways:	SNARE interactions in vesicular transport



[View online »](#)