

## Product datasheet for **TP725899**

### **LRRC59 Human Recombinant Protein**

#### **Product data:**

<b>Product Type:</b>	Recombinant Proteins
<b>Description:</b>	Recombinant Human Leucine-rich repeat-containing protein 59/LRRC59(C-His)
<b>Species:</b>	Human
<b>Expression cDNA Clone or AA Sequence:</b>	T2-S244
<b>Tag:</b>	C-His
<b>Buffer:</b>	Supplied as a 0.2 um filtered solution of 500mM PBS, 300mM NaCl, 10%Glycerol, pH7.4
<b>Note:</b>	MADE TO ORDER   Recombinant Human Leucine-rich repeat-containing protein 59/LRRC59 is produced by E.coli expression system and the target gene encoding T2-S244 is expressed with a His tag at the C-terminus.
<b>Storage:</b>	Store at $\leq -70^{\circ}\text{C}$ , stable for 6 months after receipt.Store at $\leq -70^{\circ}\text{C}$ , stable for 3 months under sterile conditions after opening.Please minimize freeze-thaw cycles.
<b>Stability:</b>	12 months from date of despatch
<b>Locus ID:</b>	55379
<b>UniProt ID:</b>	<a href="#">Q96AG4</a>
<b>Synonyms:</b>	Ribosome-binding protein p34; p34; Leucine-rich repeat-containing protein 59, N-terminally processed
<b>Protein Families:</b>	Transmembrane



[View online »](#)