

Product datasheet for **TP725903**

RFC4 Human Recombinant Protein

Product data:

| | |
|---------------------------------------|---|
| Product Type: | Recombinant Proteins |
| Description: | Recombinant Human Replication factor C subunit 4/RFC4(C-His) |
| Species: | Human |
| Expression cDNA Clone or AA Sequence: | Q2-C363 |
| Tag: | C-His |
| Buffer: | Supplied as a 0.2 um filtered solution of 500mM PBS, 300mM NaCl, 10%Glycerol, pH7.4 |
| Note: | MADE TO ORDER Recombinant Human Replication factor C subunit 4/RFC4 is produced by E.coli expression system and the target gene encoding Q2-C363 is expressed with a His tag at the C-terminus. |
| Stability: | 12 months from date of despatch |
| Locus ID: | 5984 |
| UniProt ID: | P35249 |



[View online »](#)