

Product datasheet for **TP725910**

TADA3L (TADA3) Human Recombinant Protein

Product data:

Product Type:	Recombinant Proteins
Description:	Recombinant Human Transcriptional adapter 3/TADA3(C-His)
Species:	Human
Expression cDNA Clone or AA Sequence:	S2-R369
Tag:	C-His
Buffer:	Supplied as a 0.2 um filtered solution of 500mM PBS, 300mM NaCl, 10%Glycerol, pH7.4
Note:	MADE TO ORDER Recombinant Human Transcriptional adapter 3/TADA3 is produced by E.coli expression system and the target gene encoding S2-R369 is expressed with a His tag at the C-terminus.
Stability:	12 months from date of despatch
Locus ID:	10474
UniProt ID:	<u>O75528</u>



[View online »](#)