

## Product datasheet for **TP725911**

### SMNDC1 Human Recombinant Protein

#### Product data:

Product Type:	Recombinant Proteins
Description:	Recombinant Human SMNDC1 Protein(C-His)
Species:	Human
Expression cDNA Clone or AA Sequence:	S2-Q238
Tag:	C-His
Buffer:	Supplied as a 0.2 um filtered solution of 500mM PBS, 300mM NaCl, 10%Glycerol, pH7.4
Note:	MADE TO ORDER   Recombinant Human SMNDC1 Protein is produced by E.coli expression system and the target gene encoding S2-Q238 is expressed with a His tag at the C-terminus.
Storage:	Store at $\leq -70^{\circ}\text{C}$ , stable for 6 months after receipt.Store at $\leq -70^{\circ}\text{C}$ , stable for 3 months under sterile conditions after opening.Please minimize freeze-thaw cycles.
Stability:	12 months from date of despatch
Locus ID:	10285
UniProt ID:	<a href="#">O75940</a>
Synonyms:	Survival of motor neuron-related-splicing factor 30; 30 kDa splicing factor SMNrp; SMN-related protein; Survival motor neuron domain-containing protein 1
Protein Families:	Stem cell - Pluripotency
Protein Pathways:	Spliceosome



[View online »](#)